

# Hem-o-lok® Security You Can Feel

Now available in XL

**WECK®**  
TELEFLEX MEDICAL

Hem-o-lok® provides secure, rapid ligation for open or laparoscopic surgery.

**TACTILE FEEDBACK**

Confirm closure of the Hem-o-lok® clip simply by the tactile feel.

**FLEXIBLE HINGE**

The Hem-o-lok® clip flexes in the applier, preventing clip fallout, saving time and reducing costs.

**POLYMER CLIPS**

Hem-o-lok® clips are made from a permanent, non-absorbable polymer that is inert, non-toxic and radiolucent.



**Teleflex®**  
MEDICAL

**CLIPS**

**OPEN APPLIERS/  
REMOVERS**

**ENDOSCOPIC APPLIERS/  
REMOVERS**

**AUTOMATIC  
APPLIERS**



**M**

544220  
Hem-o-lok® Medium  
Polymer Clips  
6 clips/cart, 14 cartridges per box, sterile

544113 8" Curved  
544114 11" Right-Angle (70°)  
544115 11" Curved  
544123 11" Open Remover  
(for M, ML, L)



**ML**

544230  
Hem-o-lok® Medium-Large  
Polymer Clips  
6 clips/cart, 14 cartridges per box, sterile

544170 8" Curved  
544171 11" Curved  
544172 11" Right-Angle (70°)  
544123 11" Open Remover  
(for M, ML, L)



**L**

544240  
Hem-o-lok® Large  
Polymer Clips  
6 clips/cart, 14 cartridges per box, sterile

544179 11" Right-Angle (70°)  
544180 8" Curved  
544181 11" Curved  
544123 11" Open Remover  
(for M, ML, L)



**XL**

544250  
Hem-o-lok® Extra-Large  
Polymer Clips  
6 clips/cart, 14 cartridges per box, sterile

544191 11" Curved  
544192 11" Right-Angle (70°)  
544124 11" Open Remover  
(for L, XL)

544965 Endo5® ML  
32.5cm, 5mm diameter  
544121 32cm Endo Remover  
(for ML, L)

543965 Auto Endo5®  
15 clips per applicator,  
sterile, disposable  
for endo surgery

544995 Endo10™ L  
32cm curved, 10mm diameter  
544121 32cm Endo Remover  
(for ML, L)


543181 Auto Open  
15 clips per applicator,  
sterile, disposable  
for open surgery

544990 Endo10™ XL  
32cm curved, 10mm diameter  
544130 32cm Endo Remover  
(for L, XL)

**ACCESSORIES**

544500 Sterility Storage Tray  
Storage capacity:  
3 Endoscopic instruments

**TARGET VESSEL / TISSUE BUNDLE SIZE RANGE\***



\* Hem-o-lok® ligating clips are generally suited to vessel and tissue structure sizes indicated above. Surgeon judgment should dictate clip selection for specific applications. Read instructions for use.

**SURGICAL SPECIALTIES\*\***

- General
- Colorectal
- Urology
- Gynecology

\*\*Not intended for contraceptive tubal occlusion. Surgeon judgment should dictate clip selection for specific applications.



Hem-o-lok® XL clips used during a hysterectomy.

**ORDERING OPTIONS**

Order direct from your Teleflex Medical sales professional or contact Customer Service:

Phone: 1.800.234.9325  
Fax: 1.800.932.5329  
www.teleflex.com

**INTERNATIONAL**

Phone: +919.361.3961  
Email: international@tfxsurgical.com

**TELEFLEX MEDICAL**

Teleflex Medical is a division of Teleflex Incorporated (NYSE:TFX), a diversified global organization with annual revenues of approximately \$2 billion. Teleflex designs, manufactures and distributes quality-engineered products and services for the marine, industrial, automotive, aerospace, power generation and medical markets worldwide. The company employs more than 18,000 people worldwide who focus on providing innovative solutions for customers.