


Surgical Area: Scalp

DuraHook and DermaHook Neuro Elastic Retractors and Raney Clips



Actual Size



Raney Clips

Catalog No. 658505 (clips) 10/pouch, 20 pouches/box
No. 658500 (applicators)

Raney clips are used for fast, temporary control of scalp bleeding in surgery. Weck Disposable Raney Clips effect gentle pressure with no flanges or protrusions to hamper surgical technique, and without the bulk of traditional hemostats. Weck Raney Clips may be applied with standard Raney applying forceps and are conveniently disposable.

**WECK
CLOSURE
SYSTEMS™**

One Weck Drive
Research Triangle Park, NC 27709
Phone: 800 234-9325 Fax: 800 932-5329
www.weckclosuresystems.com

World's leader in manufacturing
and development of
advanced ligating systems.