



NextStitch™
Half the Needles. Half the Time.

VALVE SUTURE ORGANIZERS



Convenient suture guides for tangle-free organization of multiple sutures used during valve implantation.

Gabbay-Frater® Suture Organizer The original valve suture organizer



4 packs per box
3 guides per pack
Ordering Information: GFS-10

4 packs per box
1 guide per pack
Ordering Information: GFS-1

Genzyme Valve Suture Organizer Improved design for smaller incisions



4 packs per box
3 guides per pack
Ordering Information: 89-9575

NEXTSTITCH™ NEEDLES



K Point Needles



Conical shaped needle
with single-cutting edge










Taper Point Needles



Round shaped needle
with tapered point

Description

Needle Length
Curvature

 KC-2	 C-2	25.9 mm 1/2 circle
 KT-5	 T-5V	19.8 mm 1/2 circle
 KT-2	—	17.8 mm 3/8 circle
 KT-3	 T-3	16.2 mm 1/2 circle

Needle illustrations are shown actual size.
NextStitch™ is a trademark of Genzyme Corporation.
Gabbay-Frater® is a registered trademark of
Genzyme Corporation.

For additional information, contact your Genzyme
representative or call 800-367-7874.

genzyme
BIOSURGERY

A division of
GENZYME CORPORATION
600 Airport Road
Fall River, MA 02720 USA
800-367-7874
508-677-6600
www.genzymebiosurgery.com

© 2001 Genzyme Corporation.
All rights reserved. Printed in USA.
Code No. 25-1380 2001-01