

2-Crank Gatch Bed with Hi-Lo Tilting (KA-462A, KA-463A)



Photograph is KA-462A

● The Hi-Lo Indicator



- The hi-lo indicator is positioned for easy reading during bed height adjustment. The Trendelenburg and Reverse Trendelenburg positions can be ascertained at a glance, ensuring that the adjustment is appropriate for treatment requirements.




- The Hi-Lo Indicator is located at both head and foot board. This makes it easy to check while making adjustment with handles.



- Crank operated back raise and knee raise adjustment
- Trendelenburg and Reverse Trendelenburg adjustment (0° - 6°)
- Hi-Lo Indicator
- 100 mm dia. Paramount original caster with brakes on opposite corners
- Pressed steel bottom

Width (mm)	
830	910
KA-462A	KA-463A



I.V.Stand (PC-5003)



- Stainless and iron main pipe set
- Dimensions (mm):
(437W x 460D x 1290-2075 H)
- 40mm dia. caster

Bedside Cabinet (PF-1100, PF-1200)



PF-1100



PF-1200

- Steel body and door
- Plastic top surface, drawer and pull-out table
- Towel rails
- Plastic shelves
- 40mm dia.caster

Type	Dimensions (mm)		
	W	D	H
PF-1100	454	469	930
PF-1200	454	469	780

PT. PARAMOUNT BED INDONESIA

MM2100 INDUSTRIAL TOWN, BLOCK M-1-1
EXPORT PROCESSING ZONE, CIKARANG BARAT
BEKASI 17520, INDONESIA
TEL: (021) 8980715 FAX: (021) 8980716
E-mail : marketing@paramount.co.id
Website : www.paramount.co.id



Reg. No. : 0853/05
Reg. No. : 04 100 20050482



PARAMOUNT BED CO.,LTD.

14-5, 2 CHOME HIGASHISUNA, KOTO-KU, TOKYO
136-8670 JAPAN

OVERSEAS TRADE DIVISION

TEL: (03) 3648-2961 FAX: (03)3648-2420

URL <http://www.paramount.co.jp>
Email : overseas@paramount.co.jp

(2005.08)

Bassinet Cart (PB-1100, PB-1200)



Photograph is PB-1100

- Clear plastic bassinet
Dimensions (mm): (755 L x 438 W x 250 H)
- Nameplate holder
- 75 mm dia. caster with brakes on opposite corners
- Polyurethane foam mattress
Dimensions (mm): (650 L x 350 W x 50 TH)
- Wire basket with powder coating finishing



Colour	
Pink	Blue
PB-1100	PB-1200



Pediatric Bed (PB-2100, PB-2200)



Photograph is PB-2100

- Crank-operated back raise adjustment
- Steel sliding side rails
- Pressed steel bottom
- 100 mm dia. Paramount original caster with brakes on opposite corners
- Polyurethane foam mattress
Dimensions (mm) : (1660 L x 830 W x 100 TH)
- I.V.Hanger Rod (Option) can be attached



Back Raise
0° - 75°

Colour	
Pink	Blue
PB-2100	PB-2200



Psychiatric Bed (PB-3101, PB-3102, PB-3103, PB-3104)



Photograph is PB-3102

- 7 pressed steel bottoms are welded to 3 long bars
- Each bottoms have holes for restraint belt
- Bed frame is supported by four legs and additional pipe leg in the middle
- Head and foot board are fixed with screws
(Board type can be chosen)

Board Type			
Steel Pipe		Plastic	
Without Corner Bumper	With Corner Bumper	Green	Beige
PB-3101	PB-3102	PB-3103	PB-3104

